



“♪ CHIRP, CHIRP ♪”

A Little Birdie Told Me About Developments in ABR & ASSR Testing;

Using the CE Chirp Stimulus for Faster Test Times

*Jill Craig, MA
2010 EHDI Meeting-Chicago
March 1, 2010*



Interacoustics[®]

leading diagnostic solutions


In a Perfect World



Interacoustics[®]

leading diagnostic solutions

In the Real World



Interacoustics[®]

leading diagnostic solutions

Challenges of Diagnosing Hearing Loss

- Can take multiple visits to finalize results
- Each appointment can take sometimes hours to complete
- Sedation or no sedation?
- Objective vs Subjective?
- “Peak Picking” Experience



Interacoustics[®]

leading diagnostic solutions

The Auditory Brainstem Response (ABR)


- An evoked potential is an electrical response induced by sound; it arises from the structures within the ear, nerve, and brain at some distance from skin electrodes.
- It is measured in the time domain
- Relies on synchronous nerve firing



Interacoustics[®]

leading diagnostic solutions

ABR...the good



- Historical, Gold Standard
- Normative Data/Research Supported
- Differential diagnosis ANSD
- Comfortable!



Interacoustics[®]

leading diagnostic solutions

ABR.....the challenges

- Can Be Time Consuming
- Peak Picking Experience
- Intensity limitations
- Single Intensity and Frequency



Interacoustics[®]

leading diagnostic solutions

The Auditory Steady State Response

- Similar to ABR, but EEG activity is analyzed in frequency domain
- Stimulus is modulated pure tone
- EEG activity modulating at same frequency as the stimulus is representative of a response
- The presence or absence of a response is determined by a statistical computer algorithm



Interacoustics[®]

leading diagnostic solutions

ASSR History

- Late 1970's and Early 1980's-Initial Research
 - Galambos
- 1991-First clinically available system
- Today-Several Clinical Systems
 - Biologic, GN Otometrics, GSI, IHS, Interacoustics



Interacoustics[®]

leading diagnostic solutions

ASSR.....the good

- More Frequency Specific
- No peak picking—Objective Detection
- Multiple and Simultaneous Frequency Testing
- Can test at higher intensities than ABR



Interacoustics[®]

leading diagnostic solutions

ASSR.....the challenges

- Limited normative/research data as compared to ABR
- Small amplitude response; maybe more sensitive to noise
- Exact Neural Generators not defined
- CHANGE!!



Interacoustics[®]

leading diagnostic solutions

Testing Considerations

- Stimulus Methods
- Detection Methods



Interacoustics[®]

leading diagnostic solutions

ABR Stimulus

- Click
 - Abrupt and rapid onset
 - Broad spectrum (theoretically stims the entire basilar membrane) *NOT FREQUENCY SPECIFIC!!!*
 - Test 2000-4000Hz frequency range!
 - The greater number of neurons that fire results in a larger response amplitude.
 - *Need Good Neural Synchrony*



Interacoustics[®]

leading diagnostic solutions

ABR Stimulus

- **Tone Burst**
 - Provides more frequency specific information
 - High correlation between behavioral and TB responses
 - Can diagnosis low and high frequency HL
 - 500Hz can be difficult!
 - Repeatability is difficult because there is less synchronous activity at that region on the cochlea
 - Requires longer window
 - Response can be 4-8ms later than a click




Interacoustics[®]

leading diagnostic solutions

New Stimulus for ABR & ASSR

Travel time in the cochlea is different for different frequencies



Interacoustics[®]

leading diagnostic solutions

Output compensation for traveling time


- If the narrow band activity is recorded, compensation for the traveling time can be obtained by time-shifting the narrow band activity
- After summation the so-called *Stacked ABR* is obtained
- The Stacked ABR is significantly larger than the normal ABR
- The procedure to obtain the Stacked ABR is complex and needs considerable recording time – and has therefore not become clinically popular.



Interacoustics[®]

leading diagnostic solutions

Stacked ABR – effective but impractical



M. Don – House Ear Institute, 2002



Interacoustics[®]

leading diagnostic solutions

Input compensation for traveling time

- Another way to compensate for the traveling time is to time-shift the different frequency components of the stimulus
- This is done by allowing the low-frequencies to appear before the high-frequencies
- Such a click with re-shuffled frequency components is called a *Chirp*.
- A *chirp* stimulus is particularly effective at lower stimulation levels , where response amplitude almost doubles compared to an equal bandwidth click stimulus.*

* Elberling and Don 2008, J. Acoust. Soc. Am.




Interacoustics[®]

leading diagnostic solutions

Instead of Synchronizing the Response....Synchronize the Stimulus!


Lower frequencies are sent a bit earlier into the cochlea
than the higher frequencies



Interacoustics[®]

leading diagnostic solutions

Frequency Specific Narrow Band CE-Chirps



Interacoustics[®]

leading diagnostic solutions

Understanding the CE Chirp

Animation of the CE Chirp



Interacoustics[®]

leading diagnostic solutions

Supporting literature

Study #1 on adults

- Testing:
 - **Delay Comp. Click** vs. **Click** (= full frequency range)
 - 49 normal-hearing younger adults







Interacoustics[®]

leading diagnostic solutions

Supporting literature

Grand Average ASSR temporal waveforms



Interacoustics[®]


leading diagnostic solutions

Supporting literature

Conclusion study #1

30 dBnHL

	Detection Rate	Detection Time
Click	83.3 %	72 s
Delay Comp. Click	97.7 %	30 s



Interacoustics[®]

leading diagnostic solutions

Supporting literature

Study #2 on newborns


- **Screening:**
 - **Click** (40 dBnHL) and **Chirp** (35 dBnHL)
 - Several studies – incl two groups of newborns (each of about N = 1,800)



Interacoustics[®]

leading diagnostic solutions

Supporting literature




- **Click - 40 dBnHL**
 - maximum test time: 120 s
 - detection criterion: 0.1 %
 - number of ears: 1744
 - detection rate: **95.4 %**
 - detection time: **42 s** (median)
47 s (mean)



Interacoustics[®]

leading diagnostic solutions

Supporting literature



- **Chirp - 35 dBnHL**
 - maximum test time: 180 s
 - detection criterion: 0.1 %
 - number of ears: 1833
 - detection rate: **96.3 %**
 - detection time: **28 s** (median)
38 s (mean)



Interacoustics[®]

leading diagnostic solutions

Supporting literature

Conclusion study #2

	Detection Rate	Detection Time
Click <u>40dBnHL</u>	95.4 %	42 s
Delay Comp. Click <u>35dBnHL</u>	96.3 %	28 s



Interacoustics[®]

leading diagnostic solutions

Documenting Efficiency of CE-Chirp[®] family

The efficiency of these new stimuli are documented in a series of publications:

1. Don, M., Elberling, C., and Maloff, E. (2009). "Input and output compensation for the cochlear traveling wave delay in wide-band ABR recordings: Implications for small acoustic tumor detection," *J. Am. Acad. Audiol.* 20, (2).
2. Elberling, C., and Don, M. (2008). "Auditory brainstem responses to a chirp stimulus designed from derived-band latencies in normal-hearing subjects," *J. Acoust. Soc. Am.* 124, 3022-3037.
3. Elberling, C, Don, M, Cebulla, M, & Stürzebecher, E. (2007). Auditory steady-state responses to chirp stimuli based on cochlear traveling wave delay. *J. Acoust. Soc. Am.* 122, 2772-2785
4. Stürzebecher, E, Cebulla, M, Elberling. C, & Berger, T. (2006). New efficient stimuli for evoking frequency-specific auditory steady-state responses. *J. Am. Acad. Audiol.* 17, 448-461.
5. Cebulla, M, Stürzebecher, E, Elberling, C, & Berger, T. (2006). New click-like stimuli for newborn hearing screening. *J. Am. Acad. Audiol.* 18, 725-738



Interacoustics[®]

leading diagnostic solutions

Implementation of ASSR into the Clinic


ASSR is an additional tool in the “Cross Check Principle”

*“In summary, we believe that the unique limitations of conventional behavioral audiometry dictate the need for a ‘test battery’ approach. The key concept governing our assessment strategy is the **cross-check principle**. The basic operation of this **principle** is that no result be accepted until it is confirmed by an independent measure. . . . We believe that the application of the cross-check **principle** to our clinical population has had an appreciable effect on the accuracy with which we can identify and quantify hearing loss during the critical years for language-learning.” (Jerger and Hayes, 1976, p. 65)*



Interacoustics[®]

leading diagnostic solutions



Interacoustics®

leading diagnostic solutions


Data Samples



Interacoustics[®]

leading diagnostic solutions

The ASSR Screen



To optimize session strategy decisions as test progresses, the response confidence is tracked over time for each test frequency.



Interacoustics[®]

leading diagnostic solutions

ASSR Audiogram

Stimulus

Right

Freq	Running	Waiting
500Hz		50dB
1kHz		50dB
2kHz		50dB
4kHz		50dB
WN		

Left

Freq	Running	Waiting
500Hz		50dB
1kHz		50dB
2kHz		50dB
4kHz		50dB
WN		

Total Session Status







Time elapsed 0:05:36

Stimulus rate

Child (30Hz)

Start

Pause



dB	500Hz	1kHz	2kHz	4kHz	dB	500Hz	1kHz	2kHz	4kHz
10	100%	56%	100%	100%	10	26%	100%	100%	100%
20	100%	100%	100%	100%	20	100%	100%	100%	100%
30	100%	100%	100%	100%	30	100%	100%	100%	100%

ASSR **Audiogram**

Stimulus

Right

Freq	Running	Waiting
500Hz		45dB
1kHz		45dB
2kHz		45dB
4kHz		45dB
WN		

Left

Freq	Running	Waiting
500Hz		45dB
1kHz		45dB
2kHz		45dB
4kHz		45dB
WN		

Total Session Status


Time elapsed 0:08:22

Stimulus rate

Adult Awake (40Hz)

Start

Pause



	500Hz	1kHz	2kHz	4kHz		500Hz	1kHz	2kHz	4kHz
20	100% 16nV	100% 12nV	100% 13nV	40% 8nV	20	100% 18nV	100% 26nV	100% 18nV	100% 13nV
30				100% 18nV	30				
40	100% 19nV	100% 15nV	100% 19nV	100% 16nV	40	100% 23nV	100% 32nV	100% 27nV	100% 32nV

They can't all be good!

- Failed UNHS (OAE and AABR)
- Full Term
- No Family History
- 5 months old at time of testing
- Normal Tymps/Absent Reflexes



Interacoustics[®]

leading diagnostic solutions

OAE Results

f2	DP	L1	L2	DP level	Noise level	S/N level
500 Hz	318 Hz	65 dB	55 dB	16.5 dB	13.0 dB	3.5 dB
1000 Hz	638 Hz	65 dB	55 dB	9.9 dB	7.2 dB	2.7 dB
2000 Hz	1278 Hz	65 dB	55 dB	9.1 dB	2.5 dB	6.6 dB
4000 Hz	2556 Hz	65 dB	55 dB	-3.2 dB	-9.4 dB	6.2 dB
6000 Hz	3836 Hz	65 dB	55 dB	-3.2 dB	-12.4 dB	9.2 dB
8000 Hz	5114 Hz	65 dB	55 dB	-12.0 dB	-9.5 dB	-2.5 dB
S/N stop criteria		Rejection level		Stimulus tolerance		
7 dB		20 dB		± 3 dB		

Right Ear

f2	DP	L1	L2	DP level	Noise level	S/N level
500 Hz	318 Hz	65 dB	55 dB	19.5 dB	14.6 dB	4.8 dB
1000 Hz	638 Hz	65 dB	55 dB	4.3 dB	-1.2 dB	5.5 dB
2000 Hz	1278 Hz	65 dB	55 dB	-8.4 dB	-5.5 dB	-1.9 dB
4000 Hz	2556 Hz	65 dB	55 dB	2.7 dB	5.0 dB	-2.2 dB
6000 Hz	3836 Hz	65 dB	55 dB	6.9 dB	5.7 dB	1.3 dB
S/N stop criteria		Rejection level		Stimulus tolerance		
7 dB		20 dB		± 3 dB		

Left Ear



Interacoustics[®]

leading diagnostic solutions

ASSR Audiogram

Stimulus Right

Freq	Running	Waiting
500Hz		45dB
1kHz		45dB
2kHz		45dB
4kHz		45dB
WN		

Left

Freq	Running	Waiting
500Hz		45dB
1kHz		45dB
2kHz		45dB
4kHz		45dB
WN		


Total Session Status
Time elapsed 0:19:52

Stimulus rate

Child (90Hz)

Start

Pause



dB	500Hz	1kHz	2kHz	4kHz	dB	500Hz	1kHz	2kHz	4kHz
20		15%			20	38%	100%		
25					25	74%			
30	51%	100%	20%	27%	30	100%	100%	67%	17%
40	32%		100%	29%	40			100%	16%
50	100%	100%	100%	77%	50	100%	100%	100%	100%
60				100%	60				

Generates an Audiogram

ASSR Audiogram

Stimulus
Right
Freq Running Waiting
500Hz
1kHz
2kHz
4kHz
WN

Left
Freq Running Waiting
500Hz
1kHz
2kHz
4kHz
WN

Total Session Status
Time elapsed 0:19:52


Stimulus rate
Child (90Hz)

Selected Correction Factor:
Standard Correction v. 1.00 (Prelim)

Start
Pause




A = ASSR result
o = Estimated Threshold Right (Default)
○ = Estimated Threshold Right (Default)



Estimated Audiogram

A = ASSR result
x = Estimated Threshold Left (Final)
X = Estimated Threshold Left (Default)



Interacoustics[®]

leading diagnostic solutions

ASSR AC Thresholds

- Many studies published and in process
 - Recently, Van Maanen & Stapells, 2009 found:

	250Hz	500Hz	1000Hz	2000Hz	4000Hz
Infants and Children		50	45	40	40

Must find this level or better to be called normal!!

MORE RESEARCH IS NEEDED!!

There is NOT sufficient evidence to support the use of ASSR only!



Interacoustics[®]

leading diagnostic solutions

ASSR BC Thresholds

- Few studies published

	500Hz	1000Hz	2000Hz	4000Hz
Preterm Infants	≤30	≤30	≤50	≤50
Post-term infants (0-11 months)	≤30	≤20	≤40	≤30
Infants (12-24 months)	≤40	≤20	≤40	≤30
Adults	≤50	≤40	≤30	≤30

Source: Recommendations from: Small, S.A. & Stapells, D.R. (2005). Multiple auditory steady-state response thresholds to bone-conduction stimuli in adults with normal hearing *Journal of the American Academy of Audiology*, 16(3): 172-183; Small, S.A & Stapells, D.R. (2006). Multiple auditory steady-state response thresholds to bone-conduction stimuli in young infants with normal hearing. *Ear and Hearing*, 27: 219-228; Small, S.A.& Stapells, D.R. (2008). Maturation of bone-conduction multiple auditory steady-state responses. *International Journal of Audiology*, 47: 476-488 ,



Interacoustics[®]

leading diagnostic solutions

Correction Factors Are Not All The Same!

System setup

Auto Tests | General Setup | Report Templates | **Correction Factors**

Predefined correction templates :

Selected correction: [Dropdown] Name of correction: [Text] [New...] [Save] [Delete...] [Cancel]

ASSR to Audiogram Correction Factors

	0/5	10/15	20/25	30/35	40/45	50/55	60/65	70/75	80/85	90/95	100	dB ASSR
500Hz	15	15	15	10	10	15	5	5	5	0	0	
1kHz	10	10	10	10	5	5	5	5	0	0	0	
2kHz	5	5	5	5	5	5	0	0	0	0	0	
4kHz	15	15	15	10	10	10	5	5	5	0	0	

[OK] [Cancel]



Interacoustics[®]

leading diagnostic solutions

Other Considerations for ASSR

- Future Research
 - Hearing Aid Fittings
 - Hearing Screening
 - Assessing Suprathreshold Hearing
 - ANSD



Interacoustics[®]

leading diagnostic solutions

One Step Further.....


- If we can generate larger response amplitudes for ASSR.....

.....how about for ABR???



Interacoustics[®]

leading diagnostic solutions



Take Home Messages

- Use of the CE Chirp is well documented in the literature
- Use of the CE Chirp can increase response amplitudes resulting in a shorter test time
- ASSR can be a part of the cross check principle



Interacoustics[®]

leading diagnostic solutions

THE END!



jgc@interacoustics-us.com



Interacoustics[®]

leading diagnostic solutions